



2335 208th Street • Torrance, CA 90501 • Tel: 310-328-1306 • Fax: 310-328-0869
<http://www.rock-ola.com>

Technical Advisory


Date: 5/2/02	Product: All SyberSonic Jukeboxes	Document # 5050
---------------------	--	------------------------

A new version 5.1d software release is being phased into production. This version replaces the 5.1.

FEATURE:

Auto Test 150 has an added timed feature that allows the user to enter the number of hours to test the jukebox.

The “TEST HOURS” can be chosen from 00 to 99. 00 is continuous test. At the end of the test, the computer will shut the jukebox down until the user presses a button on the remote or volume control. The computer will then power up and display “TEST OK” or “ERRORS EXIST”. The user can check Auto Test 151 View Report for details. If the Auto Test is being run on a C21 DigiTrack Mechanism and a fatal error occurs, the test will stop regardless of the time left.



Further information may be obtained by calling Rock-Ola Customer Service at 310-328-1306 or e-mail: service@rock-ola.com.